

# The GreenMAMBA Retractors you have always been looking for.

## Technical Information:

- 1 deployment in 4 turns
- 2 validated disinfection & sterilisation process
- 3 black shaft to prevent flares
- 4 8 and 12cm wide
- 5 8cm, articulating
- 6 5mm diameter, 33cm length



<b>REF</b>	<b>GMR12</b>	12 cm wide
<b>REF</b>	<b>GMR8B</b>	8 cm wide, articulating
<b>REF</b>	<b>GMR8</b>	8 cm wide
<b>REF</b>	<b>FCXL</b>	Flush cylinder for cleaning of the GMR12
<b>REF</b>	<b>FCL</b>	Flush cylinder long for cleaning of the GMR8(B)