


# The 1st choice in CalibrationTUBES


## Calibration



- non kinking, silicon tube
- soft tip
- cm marking
- clip to close tube
- hole in tip to avoid vacuum



TUBE-FL\_p/Rev.01



REF	<b>TUBE32W</b>	32 Fr/Ch = 10,7 mm
REF	<b>TUBE34R</b>	34 Fr/Ch = 11,3 mm
REF	<b>TUBE36B</b>	36 Fr/Ch = 12,0 mm
REF	<b>TUBE38G</b>	38 Fr/Ch = 12,7 mm
REF	<b>TUBE40Y</b>	40 Fr/Ch = 13,3 mm



www.bariatric-solutions.com  
info@bariatric-solutions.com  
Phone: +41 52 741 65 00

